

## Master - Belt Grinder

For coarse grinding of metallographic samples - used with coarse grit endless emery belt 80, 120, 220 etc., after cutting or moulding process, Removes burrs, slight cutting burns and other small metallurgical sample defects from previous cutting operations



- ½ HPAC motor
- 150 x 1220 mm belt with simplified belt releasing mechanism
- Rigid large flat bed for solid support

Model - heavy duty rough for shop floor and use before hardness inspection

### Lab Model

- Totally enclosed model for laboratory use with 100 x 915 mm belt

### Consumables-

- Abrasive Alumina Oxide Belt - Grit 120 pr Grit 220