

Patient Cots

Patient Cots – Semi Fowler, Full Fowler, 3 Functional Manual Cot/Motorized Cot, 5 Functional Manual/Motorized Cot available for rent or sale. Patient Cots – Semi fowler's bed position is the situation of a patient who is lying in bed in a turned backward position with the head of the bed at generally 30° to 45°. It is also called as 1 function cot/ Back rest cot.



These beds have a flexible back rest that can be moved by crank instrument from foot end called as “pulley”. Wheels mounted in the legs of this Semi Fowler Cot, make it simple to move starting with one place to another. It has breaks so the bed will not move by its own.

Features:

- Removable & interchangeable high quality beautiful ABS engineering plastic head panel and foot panel.
- Collapsible side railings. (Top bar Aluminum, Vertical bars SS, Bottom pipe epoxy)
- Epoxy coated mild steel frame work and 2 section perforated top.
- 125mm diameter noiseless castors (all with individual brakes)
- Provision for I.V. Rod on both sides of the bed.
- Freight saving knock-down construction.
- S. Tubular frames
- Height adjustable side railings
- Bows with leg attached
- Pre-treated and powder coated
- Castors and 2 sectional Mattress. – Optional