

Bomb Calorimeter

Bomb calorimeters are widely used for determining gross calorific values of liquid and solid samples. Automatically. User doesn't need to make any manual calculations, all calculations are done automatically displayed, printed and transferred to PC.



Bomb calorimeters are widely used for determining gross calorific values of liquid and solid samples. Automatically. User doesn't need to make any manual calculations, all calculations are done automatically displayed, printed and transferred to PC.

The controllers are specially designed & manufactured to achieve best possible measuring accuracy desired by the user.

Microprocessor based.

Automatic measurement and calculation of CALORIFIC VALUE/WATER EQUIVALENT .

- Big JUMBO LCD display.
- In-built printer
- 16 soft touch keypad
- Weigh tog tablet feedable.
- Water equivalent feedable.
- Full test report printout with date ,time & company name
- Temperature scanning resolution of 0.01 degree C.
- Alarm on firing & after the test completes
- Sensor open detection.
- Fuse wire open detection.
- No paper or printer door open detection.
- User interaction messages on LCD.
- Type of test selectable.
- Memory storage for data.
- Real Time Clock.

- Internal Data Logger with computer Interface.
- PC Software for Data Record.
- Soft Power On for full system